

## Product datasheet for **AM01157PU-N**

### Human IgG (Fc CH2 Domain specific) Mouse Monoclonal Antibody [Clone ID: 8A4]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (Fc CH2 Domain specific) Mouse Monoclonal Antibody [Clone ID: 8A4]
Clone Name:	8A4
Applications:	ELISA, FC
Recommended Dilution:	ELISA. Flow cytometry (1:50/1:100): use 10 µl of the suggested working dilution to label 10e6 cells in 100 µl.
Reactivity:	Human
Host:	Mouse
Immunogen:	Human IgG polyclonal
Isotype:	IgG1
Formulation:	PBS, pH 7.4, containing 0.09% sodium azide as a preservative State: Purified State: Liquid IgG fraction
Concentration:	lot specific
Purification:	Protein G affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)