

Product datasheet for **AM01158PU-N**

Human IgG4 (pFC) Mouse Monoclonal Antibody [Clone ID: HP6023 (also known as HO6023)]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG4 (pFC) Mouse Monoclonal Antibody [Clone ID: HP6023 (also known as HO6023)]
Clone Name:	HP6023 (also known as HO6023)
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot. Immunoblotting. Immunohistochemistry on Frozen Sections.
Reactivity:	Human
Host:	Mouse
Immunogen:	Pool of 5-6 individually purified Human IgG4 myeloma proteins
Isotype:	IgG3
Formulation:	Borate buffered saline State: Purified State: Liquid purified IgG fraction from Tissue Culture Supernatant Stabilizer: None Preservative: None
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Upon receipt, store undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.



[View online »](#)