

## Product datasheet for **AM01162PU-N**

### Influenza A (Nucleoprotein N2/N3) Mouse Monoclonal Antibody [Clone ID: A3]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	A3
Applications:	IF
Recommended Dilution:	<b>Immunofluorescence.</b>
Reactivity:	Influenza A Virus
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Influenza A.
Specificity:	This antibody recognizes the Nucleoprotein of Influenza A, binding strongly to types N2 and N3. This antibody shows no cross-reactivity with Influenza B or other respiratory viruses.
Formulation:	PBS State: Purified State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Influenza A Virus, Seasonal Flu



[View online »](#)