

Product datasheet for **AM01194PU-N**

Parainfluenza Virus 3 Mouse Monoclonal Antibody [Clone ID: B165]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B165
Applications:	ELISA, IF, IHC
Recommended Dilution:	Immunohistochemistry on frozen sections: 1/10 - 1/50 ELISA: 1/20 - 1/200 Immunofluorescence: 1/10 - 1/50
Reactivity:	Human Parainfluenza Viruses 3
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Specificity:	This antibody recognizes hemagglutinin from parainfluenza virus type 3.
Formulation:	PBS containing 0.09% Sodium Azide State: Purified State: Liquid purified Ig
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	HPIV3, HPIV-3, Parainfluenza Virus type 3



[View online »](#)