

Product datasheet for **AM01243PU-N**

Trichomonas vaginalis (Strain 1) Mouse Monoclonal Antibody [Clone ID: 1306/156]

Product data:

Product Type:	Primary Antibodies
Clone Name:	1306/156
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Trichomonas vaginalis
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Native
Specificity:	This antibody recognises <i>T. vaginalis</i> strain 1, an anaerobic parasitic protozoan and causative agent of the sexually transmitted disease Trichomoniasis.
Formulation:	PBS State: Purified State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)