

Product datasheet for **AM01264FC-N**

Hamster (Syrian) IgG Isotype Control

Product data:

Product Type:	Isotype Controls
Applications:	FC
Host:	Hamster
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Trinitrophenol - KLH
Specificity:	This antibody is suitable for use as a negative control for the measurement of non-specific binding of Syrian hamster monoclonal antibodies of isotype IgG to human, mouse or rat tissues.
Buffer:	PBS, pH 7.2 containing 0.09% sodium azide and 1% BSA
Concentration:	0.1 mg/ml
Purification:	Ion exchange chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light. Shelf life: one year from despatch.



[View online »](#)