

Product datasheet for **AM01264PU-N**

Hamster (Syrian) IgG Isotype Control

Product data:

Product Type:	Isotype Controls
Applications:	FC
Host:	Hamster
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Trinitrophenol - KLH
Specificity:	This antibody is suitable for use as a Negative Control for the measurement of non-specific binding of Syrian Hamster monoclonal antibodies of Isotype IgG to Human, Mouse or Rat tissues.
Buffer:	PBS, pH 7.2 containing 0.09% Sodium Azide
Concentration:	0.5 mg/ml
Purification:	Ion Exchange Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



[View online »](#)