

## Product datasheet for **AM01276PU-N**

### Mouse IgE Rat Monoclonal Antibody [Clone ID: LO-ME-2]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgE Rat Monoclonal Antibody [Clone ID: LO-ME-2]
Clone Name:	LO-ME-2
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Mouse
Host:	Rat
Immunogen:	Purified IgE from Balb/c mice.
Isotype:	IgG2a
Formulation:	PBS, pH 7.2 containing 0.09% Sodium Azide as preservative. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)