

Product datasheet for **AM01278PU-N**

IL13 Rat Monoclonal Antibody [Clone ID: JES10-35G12]

Product data:

Product Type:	Primary Antibodies
Clone Name:	JES10-35G12
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Rat
Isotype:	IgG2a
Clonality:	Monoclonal
Specificity:	This antibody recognises interleukin 13 (IL-13).
Formulation:	Borate buffered saline State: Purified State: Liquid purified Ig
Concentration:	lot specific
Conjugation:	Unconjugated
Gene Name:	interleukin 13
Database Link:	Entrez Gene 3596 Human P35225
Background:	IL-13 is derived from activated TH-2 cells. It inhibits proinflammatory cytokine release and down regulates macrophage activity in response to gamma interferon or bacterial LPS.
Synonyms:	IL-13, NC30



[View online »](#)