

## Product datasheet for **AM01351PU-N**

### Feline IgG Mouse Monoclonal Antibody [Clone ID: GPB2-2B1]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Feline IgG Mouse Monoclonal Antibody [Clone ID: GPB2-2B1]
Clone Name:	GPB2-2B1
Applications:	ELISA, WB
Recommended Dilution:	<b>ELISA:</b> This antibody can bind to virus-bound immune IgG and can therefore detect anti-FeLV, FIV, FCV, FHV, FCoV/FIPV and FPV IgG activity in ELISA. <b>Western Blot:</b> This antibody detects a band of approximately 50.0kDa on purified Feline IgG under reducing conditions.
Reactivity:	Feline
Host:	Mouse
Immunogen:	Purified Feline IgG
Isotype:	IgG1
Formulation:	PBS, pH 7.2 State: Purified State: Liquid purified IgG fraction from Tissue Culture Supernatant Preservative: 0.05% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)