

## Product datasheet for **AM02185PU-N**

### Guanylin (GUCA2A) (78-94) Mouse Monoclonal Antibody [Clone ID: MAK-L-G-11]

#### Product data:

|                        |   |
|------------------------|---|
| Product Type:          | Primary Antibodies  |
| Clone Name:            | MAK-L-G-11  |
| Applications:          | ELISA   |
| Recommended Dilution:  | ELISA.  |
| Reactivity:            | Human   |
| Host:                  | Mouse   |
| Isotype:               | IgG1  |
| Clonality:             | Monoclonal  |
| Immunogen:             | Synthetic human pro-Guanylin (aa 78-94)   |
| Specificity:           | This antibody detects pro-Guanylin (aa 78-94)   |
| Formulation:           | 0,1 M PBS pH 7,4<br>State: Purified<br>State: Lyophilized IgG fraction                          |
| Reconstitution Method: | Restore in aqua bidest to 1 mg/ml   |
| Purification:          | Protein G chromatography  |
| Conjugation:           | Unconjugated  |
| Storage:               | Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing. |
| Stability:             | Shelf life: One year from despatch.   |
| Gene Name:             | guanylate cyclase activator 2A  |
| Database Link:         | <a href="#">Entrez Gene 2980 Human Q02747</a>   |
| Synonyms:              | Guanylate cyclase activator 2A, GUCA2A, GUCA2   |



[View online »](#)