

## Product datasheet for **AM02185PU-S**

### Guanylin (GUCA2A) (78-94) Mouse Monoclonal Antibody [Clone ID: MAK-L-G-11]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	MAK-L-G-11
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Synthetic human pro-Guanylin (aa 78-94)
Specificity:	This antibody detects pro-Guanylin (aa 78-94)
Formulation:	0,1 M PBS pH 7,4 State: Purified State: Lyophilized IgG fraction
Reconstitution Method:	Restore in aqua bidest to 1 mg/ml
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	guanylate cyclase activator 2A
Database Link:	<a href="#">Entrez Gene 2980 Human Q02747</a>
Synonyms:	Guanylate cyclase activator 2A, GUCA2A, GUCA2



[View online »](#)