

Product datasheet for **AM02193PU-S**

Relaxin 1 (RLN1) Mouse Monoclonal Antibody [Clone ID: willi]

Product data:

Product Type:	Primary Antibodies
Clone Name:	willi
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human pro-Relaxin 1
Specificity:	This antibody detects pro-Relaxin-1.
Formulation:	PBS pH 7,4, 0,1% NaN ₃ State: Purified State: Lyophilized IgG fraction
Reconstitution Method:	Restore in aqua bidest to 1 mg/ml
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	relaxin 1
Database Link:	Entrez Gene 6013 Human P04808
Synonyms:	Relaxin H1, Prorelaxin H1, Pro-Relaxin 1, RLXH1, Relaxin A chain, Relaxin B chain
Note:	LocusID 6013



[View online »](#)