

## Product datasheet for **AM02194PU-N**

### RLN2 Mouse Monoclonal Antibody [Clone ID: kurd]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	kurd
Applications:	ELISA, FC
Recommended Dilution:	<b>ELISA.</b> <b>FACS.</b>
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human pro-Relaxin 1/2
Specificity:	This antibody detects pro-Relaxin-1/2.
Formulation:	PBS pH 7.4 containing 0.09% Sodium Azide as preservative State: Purified State: Lyophilized IgG fraction
Reconstitution Method:	Restore in aqua bidest to 1 mg/ml
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Database Link:	<a href="#">Entrez Gene 396891 Pig P01348</a>
Synonyms:	Relaxin A chain, Relaxin B chain
Note:	LocusID 6013/6019



[View online »](#)