

## Product datasheet for **AM02195PU-N**

### Relaxin 2 (RLN2) Mouse Monoclonal Antibody [Clone ID: fred]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	fred
Applications:	ELISA, FC
Recommended Dilution:	<b>ELISA.</b> <b>FACS.</b>
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human pro-Relaxin 2
Specificity:	This antibody detects pro-Relaxin-2. There was no cross reactivity obtained with Human Relaxin 1.
Formulation:	PBS pH 7.4 containing 0.09% Sodium Azide State: Purified State: Lyophilized IgG fraction
Reconstitution Method:	Restore in aqua bidest to 1 mg/ml.
Purification:	Protein G Chromatography
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	relaxin 2
Database Link:	<a href="#">Entrez Gene 6019 Human P04090</a>
Synonyms:	Relaxin H2, Prorelaxin H2, Pro-Relaxin 2, RLXH2, Relaxin A chain, Relaxin B chain
Note:	LocusID 6019



[View online »](#)