

## Product datasheet for **AM02203PU-S**

### Relaxin 1 (RLN1) Mouse Monoclonal Antibody [Clone ID: 6F1willi]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	6F1willi
Applications:	ELISA, FC
Recommended Dilution:	<b>ELISA.</b> <b>FACS.</b>
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human Relaxin 1
Specificity:	This antibody detects Relaxin-1. There was no cross reactivity obtained with Human Relaxin 2.
Formulation:	PBS pH 7.4 containing 0.09% Sodium Azide as preservative State: Purified State: Lyophilized IgG fraction
Reconstitution Method:	Restore in aqua bidest to 1 mg/ml
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	relaxin 1
Database Link:	<a href="#">Entrez Gene 6013 Human P04808</a>
Synonyms:	Relaxin H1, Prorelaxin H1, Pro-Relaxin 1, RLXH1, Relaxin A chain, Relaxin B chain
Note:	LocusID 6013



[View online »](#)