

Product datasheet for **AM02204PU-S**

Relaxin 1 (RLN1) Mouse Monoclonal Antibody [Clone ID: 8H11kurd]

Product data:

Product Type:	Primary Antibodies
Clone Name:	8H11kurd
Applications:	ELISA, FC
Recommended Dilution:	ELISA. FACS.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human Relaxin 1/2
Specificity:	This antibody detects Relaxin-1/2.
Formulation:	PBS pH 7,4, 0,1% NaN ₃ State: Purified State: Lyophilized IgG fraction
Reconstitution Method:	Restore in aqua bidest to 1 mg/ml.
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	relaxin 1
Database Link:	Entrez Gene 6013 Human P04808
Synonyms:	RLN1, Relaxin1, Relaxin-1, Prorelaxin H1, RLN2, Relaxin2, Relaxin-2, Prorelaxin H2
Note:	LocusID 6013/6019
Protein Families:	Secreted Protein



[View online »](#)