

Product datasheet for **AM02241PU-S**

Thyroglobulin (TG) Mouse Monoclonal Antibody [Clone ID: PB3]

Product data:

Product Type:	Primary Antibodies
Clone Name:	PB3
Applications:	ELISA, R
Recommended Dilution:	IRMA. ELISA.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human Thyroglobulin
Specificity:	This antibody detects Thyroglobulin
Formulation:	50 mM TRIS pH 7.4 State: Purified State: Lyophilized purified IgG fraction Preservative: 0.02% Sodium Azide
Reconstitution Method:	Restore in aqua bidest to 1 mg/ml.
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	thyroglobulin
Database Link:	Entrez Gene 7038 Human P01266
Synonyms:	Thyreoglobulin, TGN
Note:	LocusID 7038



[View online »](#)