

Product datasheet for **AM02263PU-S**

Vasopressin (AVP) Mouse Monoclonal Antibody [Clone ID: VAS 12-2-2]

Product data:

Product Type:	Primary Antibodies
Clone Name:	VAS 12-2-2
Applications:	ELISA
Recommended Dilution:	ELISA: 1 µg/ml.
Reactivity:	Human
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	Synthetic human (Arg8)-Vasopressin, bTG conjugated
Specificity:	This antibody detects Vasopressin
Formulation:	200 mM Na-Citrate/Tris, 500 mM NaCl, 0.02 % NaN ₃ , pH 7.5 State: Purified State: Lyophilized Ig fraction
Reconstitution Method:	Restore in aqua bidest to 1 mg/ml.
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	arginine vasopressin
Database Link:	Entrez Gene 551 Human P01185
Synonyms:	Arginine vasopressin, Antidiuretic Hormone, AVP
Note:	LocusID 551

[View online »](#)