

Product datasheet for **AM03005RP-N**

HLA Class I ABC Mouse Monoclonal Antibody [Clone ID: W6/32]

Product data:

Product Type:	Primary Antibodies
Clone Name:	W6/32
Applications:	FC
Recommended Dilution:	Flow Cytometry (1/50).
Reactivity:	Bovine, Feline, Human, Primate
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Membrane of human tonsil cells
Specificity:	The antibody W6/32 recognises virtually all nucleated human cells. It is a valuable reagent for analysing variations in HLA class I expression in different disease states e.g. liver disease, muscular dystrophy, inflammatory myopathy and other neuromuscular disorders. This antibody W6/32 is also suitable as a positive control for HLA tissue typing and crossmatching.
Formulation:	Provided in stabilizing phosphate buffered saline (PBS) solution containing 15mM sodium azide Label: PE State: Liquid purified Ig fraction Label: Conjugated with R-Phycoerythrin (PE) under optimum conditions
Concentration:	lot specific
Purification:	Size-Exclusion Chromatography
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light.
Stability:	Shelf life: one year from despatch.




[View online »](#)

Background:

HLA-class I major histocompatibility (MHC) antigens are intrinsic membrane glycoproteins expressed on nucleated cells and noncovalently associated with an invariant beta-2 microglobulin. They carry foreign determinants important for immune recognition by cytotoxic T cells, thus important for anti-viral and anti-tumour defence. Human HLA-class I antigens are represented by HLA-A, HLA-B and HLA-C molecules. MHC Class I molecules (MHC Class Ia) are expressed on the surface of all human cell types.

Synonyms:

MHC class I ABC, HLA Class 2 ABC

Product images:

Surface staining of human peripheral blood cells with anti-HLA-class I (W6/32) PE