

## Product datasheet for **AM03025FC-N**

### LAMP1 Mouse Monoclonal Antibody [Clone ID: B-T47]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B-T47
Applications:	FC
Recommended Dilution:	<b>Flow Cytometry:</b> Use 10 µl to label 10 <sup>6</sup> cells or 100 µl of whole blood.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	U937 cell line
Specificity:	This antibody recognises CD107a, a 100-120 kDa protein.
Formulation:	PBS Label: FITC State: Liquid purified Ig fraction Stabilizer: 5% BSA Preservative: 0.09% Sodium Azide Label: Fluoresceinisothiocyanat
Purification:	Ion Exchange Chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	lysosomal associated membrane protein 1
Database Link:	<a href="#">Entrez Gene 3916 Human P11279</a>
Synonyms:	LAMP-1, LAMPA, LGP120
Protein Families:	ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Lysosome



[View online »](#)