

## Product datasheet for **AM03032PU-L**

### CD24 Mouse Monoclonal Antibody [Clone ID: SN3]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	SN3
Applications:	FC, IHC
Recommended Dilution:	Flow Cytometry: 1/50 - 1/100; Use 10µl of the suggested working dilution to label 10e6 cells or 100µl whole blood. Immunohistochemistry on frozen sections. Recommended Positive Control: Tonsil.
Reactivity:	Chimpanzee
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Glycoproteins purified from human NALM-1 cell line. Spleen cells from immunised BALB/c mice were fused with cells of a mouse myeloma cell line.
Specificity:	This antibody recognises the CD24 cell surface antigen.
Formulation:	PBS, pH7.4 containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography on Protein G
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	CD24 molecule
Database Link:	<a href="#">P25063</a>
Background:	CD24 is a heavily glycosylated molecule of 35-45kD. CD24 is expressed by granulocytes and B lymphocytes, and by some activated T cells and T cell malignancies.



[View online »](#)

**Synonyms:** X62 heat stable antigen, HSA