

## Product datasheet for **AM05016PU-N**

### Mouse IgG2a Rat Monoclonal Antibody [Clone ID: B217R]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2a Rat Monoclonal Antibody [Clone ID: B217R]
Clone Name:	B217R
Applications:	ELISA
Recommended Dilution:	ELISA: 5µg/ml as capture.
Reactivity:	Mouse
Host:	Rat
Immunogen:	Purified IgG from Balb/c mice
Isotype:	IgG1
Formulation:	PBS, pH 7.4 containing 0.09% Sodium azide State: Aff - Purified State: Liquid purified Ig
Concentration:	lot specific
Purification:	Immunoaffinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)