

## Product datasheet for **AM05040PU-N**

### Polyethylene Glycol (PEG) Mouse Monoclonal Antibody [Clone ID: B139M]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B139M
Applications:	ELISA
Recommended Dilution:	<b>ELISA.</b>
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	This antibody binds to PEG-ylated compounds and macromolecules and to free PEG.
Formulation:	PBS containing 0.05% Sodium Azide as preservative. State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	HPLC Chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Note:	Do <u>Not</u> use in the presence of polymer detergents such as Tween-20.



[View online »](#)