

## Product datasheet for **AM05041BT-N**

### Polyethylene Glycol (PEG) Mouse Monoclonal Antibody [Clone ID: B141M]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B141M
Applications:	ELISA
Recommended Dilution:	Suitable for use in <b>ELISA</b> .
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	This antibody binds to PEG-ylated compounds and macromolecules and to free PEG.
Formulation:	PBS Label: Biotin State: Liquid purified Ig fraction. Preservative: 0.05% Sodium Azide
Concentration:	lot specific
Purification:	HPLC Chromatography.
Conjugation:	Biotin
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Note:	<b>Do <u>Not</u> use in the presence of polymer detergents such as Tween-20.</b>



[View online »](#)