

Product datasheet for **AM05073PU-N**

Salmonella (Serogroups A, B, C1, C2, D, E1, E3, E4, F, G1, G2) Mouse Monoclonal Antibody [Clone ID: B395M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B395M
Applications:	ELISA, IF
Recommended Dilution:	Suitable for use in ELISA and IFA.
Reactivity:	Salmonella species
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Heat-inactivated, mixture of Salmonella serogroups A, B, C1, C2, and D
Specificity:	This antibody is specific for common LPS core of all Salmonella tested. Positive for all Salmonella O-serogroups tested: A, B, C1, C2, D, E1, E3, E4, F, G1, G2. Negative against related bacteria: E.coli (mix), Klebsiella, Citrobacter, Pseudomonas, Yersinia, Shigella, Proteus and Legionella.
Formulation:	0.01M PBS, pH 7.2 containing 0,09% Sodium Azide as preservative. State: Purified State: Liquid purified Ig fraction (> 90% pure).
Concentration:	lot specific
Purification:	Protein A Chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)