

Product datasheet for **AM05119PU-N**

TSH beta (TSHB) Mouse Monoclonal Antibody [Clone ID: 057-11006]

Product data:

Product Type:	Primary Antibodies
Clone Name:	057-11006
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human pituitary gland
Specificity:	This antibody reacts to Thyroid Stimulating Hormone (TSH). Crossreacts with LH and FSH <0.10%. Does not cross react with hCG.
Formulation:	10mM Sodium Phosphate, pH 7.4 + 0.2 containing 150mM Sodium chloride and 0,09% sodium azide State: Purified State: Liquid purified Ig (0.2µm filtered)
Concentration:	lot specific
Purification:	>90% pure (SDS-PAGE). Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one week or (in aliquots) at -40°C for longer. If aliquoted for long term storage, fill volume should be equal to or greater than 50% of the nominal fill volume of the vial used. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	thyroid stimulating hormone beta
Database Link:	Entrez Gene 7252 Human P01222



[View online »](#)

- Background:** Various hormones are secreted from the anterior pituitary during development and growth. Thyroid Stimulating Hormone (TSH) is secreted from the anterior pituitary and controls thyroid structure and metabolism.
- Synonyms:** Thyrotropin subunit beta, Thyroid-stimulating hormone subunit beta, TSH, TSHB, TSH beta