

Product datasheet for **AM05226PU-N**

CDK6 Mouse Monoclonal Antibody [Clone ID: DCS-83.1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	DCS-83.1
Applications:	IHC, IP, WB
Recommended Dilution:	Western Blot. Immunoprecipitation. Immunohistochemistry.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant cdk6 protein recognizing a protein of 40kD recognizing a 40kD protein.
Formulation:	PBS containing 0.08% Sodium Azide as preservative. State: Purified State: Liquid (sterile filtered) purified IgG fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	cyclin-dependent kinase 6
Database Link:	Entrez Gene 1021 Human Q00534
Background:	It plays a crucial role in the progression of G1 to S. It can be highly expressed in some lymphoid tumor and sarcoma cell lines (ML-1 and U-2-OS resp.).
Synonyms:	PLSTIRE



[View online »](#)