

## Product datasheet for **AM05458PU-N**

### Osteocalcin (BGLAP) Mouse Monoclonal Antibody [Clone ID: 8H12F9H10]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	8H12F9H10
Applications:	ELISA
Recommended Dilution:	ELISA: This antibody may be used in combination with AM01190PU-N or AM05459PU-N in sandwich ELISA assays.
Reactivity:	Bovine, Human, Porcine, Rabbit, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant human osteocalcin GST fusion protein
Specificity:	This antibody recognizes intact human and bovine osteocalcin in a "2 site" assay and reacts with osteocalcin fragments corresponding to amino acids 1-19, 7-19 and 15-31 of the native molecule. It does not recognize fragments representing amino acids 20-43 or intact bovine osteocalcin.
Formulation:	Phosphate buffered saline containing 0.09% Sodium Azide (NaN <sub>3</sub> ) State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography on Protein G
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	bone gamma-carboxyglutamate protein
Database Link:	<a href="#">Entrez Gene 632 Human P02818</a>



[View online »](#)

**Background:** Osteocalcin (Bone Gla Protein) is a 49 amino acid single chain vitamin K dependent protein (MW 5.8 kDa), made by osteoblasts that is a major component of the non-collagenous bone matrix. Serum osteocalcin is elevated in diseases characterized by increased bone turnover such as osteoporosis, hyperparathyroidism and Pagets disease, and low in conditions associated with low bone turnover such as hypoparathyroidism and growth hormone deficiency.

**Synonyms:** Bone Gla protein, BGLAP