

Product datasheet for **AM05719PU-N**

Cytomegalovirus / CMV Late Antigen (65 kDa) Mouse Monoclonal Antibody [Clone ID: 0896]

Product data:

Product Type:	Primary Antibodies
Clone Name:	0896
Applications:	ELISA, IF, WB
Recommended Dilution:	Immunofluorescence. ELISA. Western Blot.
Reactivity:	Cytomegalovirus
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Specificity:	This antibody reacts to Cytomegalovirus 65kDa late Antigen.
Formulation:	Phosphate buffered saline containing 0.09% Sodium Azide (NaN ₃) State: Purified State: Liquid purified IgG
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	HCMV, HHV5



[View online »](#)