

## Product datasheet for **AM05722PU-N**

### **E. coli K99 Pili Mouse Monoclonal Antibody [Clone ID: 6911]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	6911
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Escherichia coli
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	This antibody reacts to <i>Escherichia Coli</i> K99 Pili.
Formulation:	Phosphate buffered saline containing 0.09% Sodium Azide (NaN <sub>3</sub> ) State: Purified State: Liquid purified IgG
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Escherichia coli, EHEC, EPEC, ETEC



[View online »](#)