

## Product datasheet for **AM05765PU-N**

### Rabies Virus Mouse Monoclonal Antibody [Clone ID: 6041]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	6041
Applications:	ELISA, IF, IHC
Recommended Dilution:	Immunofluorescence. ELISA. Immunohistochemistry on frozen sections.
Reactivity:	Rabies Lyssavirus
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Specificity:	This antibody recognises a Rabies virus surface glycoprotein.
Formulation:	Phosphate buffered saline containing 0.09% Sodium azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)