

## Product datasheet for **AM05769PU-N**

### Rotavirus (capsid protein) Mouse Monoclonal Antibody [Clone ID: 0521]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	0521
Applications:	ELISA, WB
Recommended Dilution:	Western Blot. ELISA.
Reactivity:	Rotavirus
Host:	Mouse
Isotype:	IgG2
Clonality:	Monoclonal
Specificity:	This antibody recognises the p41 capsid protien of Rotavirus from human, monkey and porcine isolates.
Formulation:	Phosphate buffered saline containing 0.09% Sodium azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)