

Product datasheet for **AM05802PU-N**

Influenza A H1N1 (Haemagglutinin) Mouse Monoclonal Antibody [Clone ID: 58AB7-19-18]

Product data:

Product Type:	Primary Antibodies
Clone Name:	58AB7-19-18
Applications:	IF
Recommended Dilution:	Immunofluorescence: 1:25 - 1:100.
Reactivity:	Influenza A Virus
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Influenza blend
Specificity:	<p>This antibody is specific for the Influenza A H1 antigen and recognizes most H1N1 strains isolated in California since 1986, including A/Taiwan/1/86 and A/USSR/90/77. Some cross reactivity has been observed with A/Texas/36/91 and the New Caledonia strain seen after 1998.</p> <p>This antibody does not recognize the following H3N2 strains: A/Wuhan/396/95, A/Johannesburg/33/94, A/Shandong/93, A/Shanghai/16/89, all CA strains isolated from 1986 to the present.</p>
Formulation:	Phosphate buffered saline containing 0.09% Sodium azide State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2 - 8 °C up to one month or (in aliquots) at -20 °C for longer. Avoid repeated freezing and thawing. Should this product contain a precipitate we recommend microcentrifugation before use.
Stability:	Shelf life: one year from despatch.
Synonyms:	Seasonal Flu H1N1



[View online »](#)