

Product datasheet for **AM06012PU-N**

IL6 Mouse Monoclonal Antibody [Clone ID: 6D9A1,6C12B1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	6D9A1,6C12B1
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA: 1/10000. Immunohistochemistry on Paraffin Sections: 1/200. Western Blot: 1/200 - 1/1000.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Ni-NTA purified recombinant Human IL-6 expressed in E. Coli strain BL21 (DE3)
Specificity:	This antibody detects Interleukin 6 (IL6) (Family: Hormone, Subfamily: not assigned-Hormone).
Formulation:	State: Ascites State: Ascitic fluid containing 0.03% Sodium Azide as preservative
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	interleukin 6
Database Link:	Entrez Gene 3569 Human P05231



[View online »](#)

Background:

Interleukin 6 (IL6) is a multifunctional 24 kDa protein originally discovered in the medium of RNA stimulated fibroblastoid cells. It is upregulated by IL1, TNF, PDGF, IFN beta, TNF alpha, NGF, IL17 and downregulated by glucocorticoids IL4, TGF beta. IL6 appears to be directly involved in the responses that occur after infection and cellular injury, and it may prove to be as important as IL1 and TNF alpha in regulating the acute phase response. IL6 has also been implicated in regulating adipose mass.

IL6 is reported to be produced by fibroblasts, activated T cells, activated monocytes or macrophages and endothelial cells. It acts upon a variety of cells including fibroblasts, myeloid progenitor cells, T cells, B cells and hepatocytes. In addition, IL6 appears to interact with IL2 in the proliferation of T lymphocytes. IL6 potentiates the proliferative effect of IL3 on multipotential hematopoietic progenitors.

Synonyms:

IL-6, Interferon beta-2, IFNB2, B-cell stimulatory factor 2, BSF-2, CDF


Product images:

Figure 1. AM06012PU-N IL6 antibody staining of Formalin-Fixed Paraffin-Embedded Human Tonsil.