

Product datasheet for **AM06072PU-N**

Human IgG Mouse Monoclonal Antibody [Clone ID: 4D2D9G8]


Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG Mouse Monoclonal Antibody [Clone ID: 4D2D9G8]
Clone Name:	4D2D9G8
Applications:	ELISA, WB
Recommended Dilution:	Western Blot: 1/500 - 1/2000. ELISA: Propose dilution 1/10000.
Reactivity:	Human
Host:	Mouse
Immunogen:	Human IgG was isolated from Human sera and purified by chromatography.
Isotype:	IgG1
Formulation:	PBS State: Purified State: Liquid purified IgG fraction Preservative: 0.03% Sodium Azide
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)

Product images:



Western blot analysis using Human IgG antibody Cat-No AM06072PU-N (Lane 1, 2) and Albumin mouse mAb (Lane 3, 4) against Human Serum (Lane 1, 3) and plasma (Lane 2, 4).