

## Product datasheet for **AM08012BT-N**

### Mouse IgD (Delta heavy chain) Rat Monoclonal Antibody [Clone ID: 11-26]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgD (Delta heavy chain) Rat Monoclonal Antibody [Clone ID: 11-26]
Clone Name:	11-26
Applications:	FC
Recommended Dilution:	Flow Cytometry: < / = 1µg/10e6 cells. (See Ref.1-5)
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: Biotin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)