

## Product datasheet for **AM08012LE-N**

### Mouse IgD (Delta heavy chain) Rat Monoclonal Antibody [Clone ID: 11-26]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgD (Delta heavy chain) Rat Monoclonal Antibody [Clone ID: 11-26]
Clone Name:	11-26
Applications:	FC, FN
Recommended Dilution:	Flow Cytometry (See Ref.1-4). Functional Assays.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a
Formulation:	PBS without preservatives. State: Low Endotoxin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. <b>DO NOT FREEZE!</b> Avoid repeated freezing and thawing.



[View online »](#)