

Product datasheet for **AM08012PU-N**

Mouse IgD (Delta heavy chain) Rat Monoclonal Antibody [Clone ID: 11-26]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgD (Delta heavy chain) Rat Monoclonal Antibody [Clone ID: 11-26]
Clone Name:	11-26
Applications:	FC
Recommended Dilution:	Flow Cytometry. (See Ref.1-5)
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a
Formulation:	100 mM Borate buffered saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)