

Product datasheet for **AM08012RP-N**

Mouse IgD (Delta heavy chain) Rat Monoclonal Antibody [Clone ID: 11-26c]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgD (Delta heavy chain) Rat Monoclonal Antibody [Clone ID: 11-26c]
Clone Name:	11-26c
Applications:	FC
Recommended Dilution:	Flow Cytometry: $< / = 0,2 \mu\text{g}/10\text{e}6 \text{ cells}$. (See Ref.1-5)
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction. Label: R-phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)