

## Product datasheet for **AM08012RP-S**

### Mouse IgD (Delta heavy chain) Rat Monoclonal Antibody [Clone ID: 11-26]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgD (Delta heavy chain) Rat Monoclonal Antibody [Clone ID: 11-26]
Clone Name:	11-26
Applications:	FC
Recommended Dilution:	Flow Cytometry: $< / = 0,2 \mu\text{g}/10\text{e}6 \text{ cells}$ . (See Ref.1-5)
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction Label: R-phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)