

Product datasheet for AM08013AP-N

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Mouse IgE (Epsilon heavy chain) Rat Monoclonal Antibody [Clone ID: 23G3]

Product data:

Product Type: Secondary Antibodies

Product Name: Mouse IgE (Epsilon heavy chain) Rat Monoclonal Antibody [Clone ID: 23G3]

Clone Name: 23G3
Applications: ELISA

Recommended Dilution: ELISA: 1/500-1/1000.

Reactivity: Mouse
Host: Rat
Isotype: IgG1

Formulation: 50mM Tris/1mM MgCl2/50% Glycerol, pH 8.0 containing 0.09% Sodium Azide as preservative.

Label: AP

State: Liquid purified lg fraction. Label: Alkaline Phosphatase

Conjugation: AF

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing.