

## Product datasheet for AM08013BT-N

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Mouse IgE (Epsilon heavy chain) Rat Monoclonal Antibody [Clone ID: 23G3]

**Product data:** 

**Product Type:** Secondary Antibodies

**Product Name:** Mouse IgE (Epsilon heavy chain) Rat Monoclonal Antibody [Clone ID: 23G3]

Clone Name: 23G3

**Applications:** ELISA, FC, IF

**Recommended Dilution:** Flow Cytometry: < / = 1  $\mu$ g/10e6 cells.

Identification and enumeration of IgE+ cells.

Immunofluorescence Microscopy: Identification and enumeration of IgE+ cells.

ELISA: 1/5000-1/10000.

Enzyme-Linked-Immunosorbent-Assay.

Reactivity: Mouse Host: Rat

Isotype:

**Formulation:** PBS containing 0.09% Sodium Azide as preservative.

Label: Biotin

lgG1

State: Liquid purified Ig fraction.

**Concentration:** lot specific

Conjugation: Biotin

**Storage:** Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

