

Product datasheet for **AM08013FC-N**

Mouse IgE (Epsilon heavy chain) Rat Monoclonal Antibody [Clone ID: 23G3]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgE (Epsilon heavy chain) Rat Monoclonal Antibody [Clone ID: 23G3]
Clone Name:	23G3
Applications:	FC, IF
Recommended Dilution:	Flow Cytometry: < / = 1 µg/10e6 cells. Identification and enumeration of IgE+ cells. Immunofluorescence Microscopy: Identification and enumeration of IgE+ cells.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: FITC State: Liquid purified Ig fraction. Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)