

Product datasheet for **AM08013LE-N**

Mouse IgE (Epsilon heavy chain) Rat Monoclonal Antibody [Clone ID: 23G3]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgE (Epsilon heavy chain) Rat Monoclonal Antibody [Clone ID: 23G3]
Clone Name:	23G3
Applications:	FC, IF
Recommended Dilution:	Functional Assays.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Formulation:	PBS without preservatives. State: Low Endotoxin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.



[View online »](#)