

## Product datasheet for **AM08013RP-N**

### Mouse IgE (Epsilon heavy chain) Rat Monoclonal Antibody [Clone ID: 23G3]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgE (Epsilon heavy chain) Rat Monoclonal Antibody [Clone ID: 23G3]
Clone Name:	23G3
Applications:	FC, IF
Recommended Dilution:	<b>Flow Cytometry:</b> < / = 0.2 µg/10e6 cells. Identification and enumeration of IgE+ cells. <b>Immunofluorescence Microscopy:</b> Identification and enumeration of IgE+ cells.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Formulation:	Label: PE State: Liquid purified Ig fraction. Label: R-phycoerythrin
Concentration:	lot specific
Purification:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent.
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)