

Product datasheet for **AM08015AP-N**

Mouse IgM (Mu chain specific) Rat Monoclonal Antibody [Clone ID: 1B4B1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgM (Mu chain specific) Rat Monoclonal Antibody [Clone ID: 1B4B1]
Clone Name:	1B4B1
Applications:	ELISA
Recommended Dilution:	ELISA: 1/1000-1/2000. Also used as a second step reagent for mouse IgM monoclonal antibodies.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Formulation:	50mM Tris/1mM MgCl ₂ /50% Glycerol, pH 8.0 containing 0.09% Sodium Azide as preservative. Label: AP State: Liquid purified Ig fraction. Label: Alkaline Phosphatase
Conjugation:	AP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)