

Product datasheet for **AM08015HR-N**

Mouse IgM (Mu chain specific) Rat Monoclonal Antibody [Clone ID: 1B4B1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgM (Mu chain specific) Rat Monoclonal Antibody [Clone ID: 1B4B1]
Clone Name:	1B4B1
Applications:	ELISA
Recommended Dilution:	ELISA: 1/2000-1/4000. Also used as a second step reagent for mouse IgM monoclonal antibodies.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Formulation:	50% PBS/50% Glycerol, pH 7.4. without preservatives. Label: HRP State: Liquid purified Ig fraction. Label: Horseradish Peroxidase
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)