

Product datasheet for **AM08017AP-N**

Mouse IgG2a (gamma 2a chain specific) Rat Monoclonal Antibody [Clone ID: SB84a]

Product data:

| | |
|-----------------------|---|
| Product Type: | Secondary Antibodies |
| Product Name: | Mouse IgG2a (gamma 2a chain specific) Rat Monoclonal Antibody [Clone ID: SB84a] |
| Clone Name: | SB84a |
| Applications: | ELISA |
| Recommended Dilution: | ELISA: 1/500-1/1000. Also used as a second step reagent for mouse IgG2a monoclonal antibodies. |
| Reactivity: | Mouse |
| Host: | Rat |
| Isotype: | IgG1 |
| Formulation: | 50mM Tris/1mM MgCl ₂ /50% Glycerol, pH 8.0 containing 0.09% Sodium Azide as preservative. Label: AP State: Liquid purified Ig fraction. Label: Alkaline Phosphatase |
| Conjugation: | AP |
| Storage: | Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. |



[View online »](#)