

Product datasheet for **AM08017BT-N**

Mouse IgG2a (gamma 2a chain specific) Rat Monoclonal Antibody [Clone ID: SB84a]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2a (gamma 2a chain specific) Rat Monoclonal Antibody [Clone ID: SB84a]
Clone Name:	SB84a
Applications:	ELISA, FC, IF
Recommended Dilution:	Flow Cytometry: Identification and enumeration of IgG2a+. Immunofluorescence Microscopy: Identification and enumeration of IgG2a+. ELISA. Also used as a second step reagent for mouse IgG2a monoclonal antibodies.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: Biotin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)